## **INTRODUCTION**

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## INTRODUCTION

The presentation of an implant card to a patient fitted with an implantable medical device is a new requirement of Regulation (EU) 745/2017. This responsibility is shared between the manufacturer and the health institution that is in charge of implanting the medical device.

The implant card aims to achieve four main objectives (items extracted from the MDR article 18 paragraph 1):

*“a) information allowing the identification of the device, including the device name, serial number, lot number, the UDI, the device model, as well as the name, address and the website of the manufacturer.*

*b) any warnings, precautions or measures to be taken by the patient or a healthcare professional with regard to reciprocal interference with reasonably foreseeable external influences, medical examinations or environmental conditions.*

*c) any information about the expected lifetime of the device and any necessary follow-up.*

*d) any other information to ensure safe use of the device by the patient, including the information in point (u) of Section 23.4 of Annex I: in the case of implantable devices, the overall qualitative and quantitative information on the materials and substances to which patients can be exposed;”*

## ROLE OF THE HEALTH INSTITUTION

There is an implant card in the packaging of each device. Health institutions have the responsibility to complete the implant card, to present it to the patient and to give him/her the following information (items extracted from the new European Regulation 2017/745 Article 18):

*“Member States shall require health institutions to make the information referred to in paragraph 1 available, by any means that allow rapid access to that information, to any patients who have been implanted with the device, together with the implant card, which shall bear their identity. “*

The patient will be able to find the information and information updates for their card on the following internet link: <http://info.xnov.com>

## HOW SHOULD THE HEALTH INSTITUTION USE THE IMPLANT CARD?

An implant card is supplied in the packaging of each implantable device.

Traceability labels are supplied in the packaging of each implantable device. The institution must complete an implant card for each procedure and affix a traceability label for each implanted device following the procedure below:

### 1- COMPLETE THE BLANK FIELDS

The card contains blank fields that must be completed by the health institution:

- 1**  Patient name or identification
- 2**  Date of implantation
- 3**  Name and address of the health institution that carried out the surgery

### 2- TICK THE TYPE OF PROSTHESIS USED

- Tick: “hip prosthesis” or “knee prosthesis” depending on the joint involved.

### 3- AFFIX THE TRACEABILITY LABELS

Apply a label from each of the implant boxes used on the card, in the places provided for this purpose.

<div style="display: flex; flex-direction: column; align-items: flex-start;"> <div style="margin-bottom: 10px;"> <b>1</b> _____</div> <div style="margin-bottom: 10px;"> <b>2</b> _____</div> <div style="margin-bottom: 10px;"> <b>3</b> _____</div> <div style="margin-bottom: 10px;">_____</div> <div style="margin-bottom: 10px;">_____</div> <div style="margin-bottom: 10px;">_____</div> </div> <div style="margin-top: 10px;"> <p><input checked="" type="checkbox"/> FR - Prothèse de hanche    <input type="checkbox"/> FR - Prothèse de genou</p> <p>EN - Hip prosthesis            EN - Knee prosthesis</p> <p>DE - Hüftprothese            DE - Knieprothese</p> <p>ES - Prótesis de cadera        ES - Prótesis de rodilla</p> <p>NL - Heuprothese            NL - Knieprothese</p> <p>IT - Protesi d'anca            IT - Protesi di ginocchio</p> <p>SK - Kyčelní protéza        SK - Kolenní protéza</p> <p>RO - Proteză de șold        RO - Proteza genunchiului</p> <p>CZ - Protéza bedrového klbu    CZ - Protéza kolena</p> </div> <div style="margin-top: 10px;"> <p> <b>X·NOV MEDICAL TECHNOLOGY</b></p> <p>Rue d'Airmont 7</p> <p>CH – 2900 PORRENTROY</p> <p>+41 32 544 65 60</p> <p> <a href="http://www.xnov.com">www.xnov.com</a></p> </div>	<div style="border: 1px dashed gray; padding: 10px; margin-bottom: 10px;">  <p><i>FR - Coller l'étiquette ici</i>  <i>EN - Stick the label here</i>  <i>DE - Hier die Etiketle aufkleben</i>  <i>ES - Pegar la etiqueta aquí</i>  <i>NL - Plak het etiket hier</i>  <i>IT - Attaccare l'etichetta qui</i>  <i>SK - Nalepte štítok sem</i>  <i>RO - Lipiți eticheta aici</i>  <i>CZ - Do tohoto místa nalepte štítek</i></p> </div> <div style="border: 1px dashed gray; padding: 10px;">  <p><i>FR - Coller l'étiquette ici</i>  <i>EN - Stick the label here</i>  <i>DE - Hier die Etiketle aufkleben</i>  <i>ES - Pegar la etiqueta aquí</i>  <i>NL - Plak het etiket hier</i>  <i>IT - Attaccare l'etichetta qui</i>  <i>SK - Nalepte štítok sem</i>  <i>RO - Lipiți eticheta aici</i>  <i>CZ - Do tohoto místa nalepte štítek</i></p> </div>
---	--

**4- GIVETHE CARD AND THE INFORMATION TO THE PATIENT**

- **Give the implant card duly completed to the patient**
- **Give the supplementary information to the patient:**

*“ b) any warnings, precautions or measures to be taken by the patient or a healthcare professional with regard to reciprocal interference with reasonably foreseeable external influences, medical examinations or environmental conditions.*

*c) any information about the expected lifetime of the device and any necessary follow-up.*

*d) any other information to ensure safe use of the device by the patient, including the information in point (u) of Section 23.4 of Annex I: in the case of implantable devices, the overall qualitative and quantitative information on the materials and substances to which patients can be exposed;”*

**SYMBOLS**

The symbols that may appear on the implant card and/or the traceability labels are defined in the table below:

	Name of the patient
	Date of the surgery
	Name and address of the health institution
	Information website for patients
	Manufacturer
	Name of the device
	Reference of the device
	Lot number
	Materials used
	Unique Device Identification
	Use-by date
	MRI-conditional compatibility

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